

General Care of Mattresses

The mattresses supplied are generally Gold Deluxe, Silver or Concave Mattresses.

These are Single Sided mattresses and must not be turned over.

The fabric cover is a two way stretch polyurethane and they have fully welded covers or in some cases a waterproof zip-on cover.

Care Instructions:

- Wipe the surface clean with a sponge and warm soapy water at up to 95°C
- Do not use abrasive cleaners, solvents or bleaches
- Do not iron

There are air holes in the ends of the mattresses. When washing mattresses, please be careful not to let water into these holes.



Misuse will void supplier's warranty on furniture